

powered by **SEGWAY®**



Transport

The reclinable backrest and the with-one-hand-removable handlebar reduce Genny 2.0 to a 67 cm cube, compact and suitable for transportation.



GM Design



Italian design
Swiss precision

GENNY®
MOBILITY



Technical Info Genny 2.0

- Range up to 15-25 km depending on driving style and road surface
- Ground clearance 9cm
- Uphill and downhill up to 18%
- Battery recharge time 6.5 hours
- Adjustable seat and backrest individually

www.gennymobility.com
info@gennymobility.com

Mobility Supplies NZ
222 Ruahine Street
Palmerston North
0800 000 605
sales@mobilitysuppliesnz.co.nz
www.mobilitysuppliesnz.co.nz

GENNY FACTORY SA
GLOBAL HEADQUARTERS
VIA ESSAGRA 7
6592 S.ANTONINO (SWITZERLAND)

Tel +41 (91) 960 70 50 | info@gennymobility.com

ENG. BRG. 030918



Rediscover your life
www.gennymobility.com



Infokey

The wireless infokey controller gives access to all the functions of Genny 2.0.

Security system

The side handles, besides giving support during use and assisting in the transfer to and from another wheelchair, trigger the mechanism of descent of the support feet in case of emergency.



Nothing is left to chance

The presence sensor under the seat allows the functioning of Genny 2.0 only if a driver is on board, avoiding accidental activation and ensuring deactivation when getting off.



White



Silver



Matte black



Glossy red



Everything under control

Control panel ON/OFF
Parking and lights (optional) control
Battery charge level indicator

GENNY[®]
MOBILITY



Parking

The parking feet allow you to place your Genny 2.0 in a safe and easy way.